

IMMUNOHISTOCHEMICAL STAINING PNA-HRP PROTOCOL



Uniformed Services University
Department of Microbiology and Immunology
4301 Jones Bridge Road
Bethesda, MD 20814

Used by the Laboratory of William C. Gause, Ph.D.

TISSUE PREPARATION FOR SECTIONING

1. Label 15mL aluminum liquid N₂ containers (Accurate Chemical Scientific, 1-800-645-6264) containing 4mL of frozen tap H₂O.
2. Sacrifice animals and drop the desired tissues in the aluminum container filled with liquid N₂. Keep groups separated. When liquid N₂ has evaporated, cap the container and store at -70°C.
3. To begin sectioning, mount the individual frozen tissue with tissue freezing medium (Scientific Products, 1-800-964-5227) and allow to freeze thoroughly.
4. Routine sections are trimmed at 20 µm and cut for slides at 5 to 8 µm and picked up on a labeled glass slide.

PNA-HRP STAIN

1. Remove frozen slides from the -70°C freezer and allow to warm to room temperature.
2. To contain the staining solutions, encircle the section including the tissue medium with a hydrophobic slide marker (PAP pen, Research Products International Corp., 1-800-323-9814).
3. Rinse 2-3X in 1X PBS.
4. Add 250µL of PNA-HRP (1:75 dilution in 1%BSA/PBS). (PNA-HRP, # L-7759, Sigma Chemical, 800.325.3010)
5. Rinse 2-3X in 1X PBS.
6. Add 250µL formamide to 4mg 3-Amino-9-Ethyl-Carbazole (AEC, Sigma: A-5754) and place in 9.75 mL sodium acetate buffer (.05M CH₃COONa-3H₂O, pH 5.0). Filter this solution (0.22µm) and add 50 µL of 3% hydrogen peroxide solution. Immediately add to the slides.
7. Developing takes 5-15 minutes and can be monitored under the microscope.
8. Rinse with 1X PBS 4 times for a minimum of 5 minutes each, longer aids in removing background.
9. To mount, apply a thin layer of crystal mount, being careful to NOT touch the tissue. Cover when drying to avoid dust.